

Anti-Ucp2 Rabbit Polyclonal Antibody For Wb, Elisa - P55851

Item Number: CM0025067

Introduction

Rabbit polyclonal antibody targeting UCP2 (Uncoupling Protein 2), validated for Western blot and ELISA, reacting with human, mouse, and rat. Research use for mitochondrial metabolism, oxidative stress, and inflammation studies.

[Learn More](#)

Parameter	Value
Product Name	UCP2 Rabbit pAb
Remarks/Alias	UCPH; BMIQ4; SLC25A8; UCP2
Species	Human
GeneID (Human)	7351
GeneID	7351
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 112-211 of human UCP2 (NP_003346.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55851
Protein Weight	33kDa
Shipping	Ice bag