

Anti-ULK1 Polyclonal Antibody For Wb, If/Icc, Elisa - 075385

Item Number: CM0028314



Introduction

ULK1 Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA in human, mouse, rat. Recognizes the 113kDa autophagy kinase ULK1. Suitable for autophagy research. UniProt 075385.

[Learn More](#)

Parameter	Details
Product Name	ULK1 Rabbit pAb
Remarks/Alias	ATG1, ATG1A, UNC51, hATG1, Unc51.1, ULK1
Species	Human
GeneID (Human)	8408
GeneID	8408
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 350-450 of human ULK1 (NP_003556.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075385
Protein Weight	113kDa
Shipping	Ice bag