

Anti-ULK4 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q96C45

Item Number: CM0025289



Introduction

Rabbit polyclonal antibody targeting human ULK4 pseudokinase, validated for WB, IF/ICC, and ELISA applications. Detects ULK4 in human, mouse, and rat. Suitable for neuroscience and hedgehog signaling research.

[Learn More](#)

Parameter	Value
Product Name	ULK4 Rabbit pAb
Remarks/Alias	FAM7C1; REC01035; ULK4
Species	Human
GeneID (human)	54986
GeneID	54986
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-290 of human ULK4 (NP_060356.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96C45
Protein Weight	142kDa
Shipping	Ice bag