


Anti-Upf3A Polyclonal Antibody For Wb, If/Icc, Elisa - Q9H1J1

Item Number: CM0017955



Introduction

Rabbit polyclonal antibody against human UPF3A, suitable for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Detects a 55 kDa band; ideal for nonsense-mediated mRNA decay and RNA transport research.

[Learn More](#)

Parameter	Value
Product Name	UPF3A Rabbit pAb
Remarks/Alias	UPF3; HUPF3A; RENT3A; UPF3A
Species	Human
GeneID (human)	65110
GeneID	65110
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 327-476 of human UPF3A (NP_075387.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H1J1
Protein Weight	55kDa
Shipping	Ice bag