

Anti-Upp1 Rabbit Monoclonal Antibody For Wb, Elisa - Q16831

Item Number: CM0002004



Introduction

Anti-UPP1 rabbit monoclonal antibody for WB and ELISA. Validated in Human, Mouse, Rat. Targets uridine phosphorylase 1, a key enzyme in pyrimidine metabolism, ideal for nucleotide metabolism research.

[Learn More](#)

Parameter	Value
Product Name	UPP1 Rabbit mAb
Remarks/Alias	UP; UPP; UPASE; UDRPASE; UPP1
Species	Human
Gene ID (Human)	7378
Gene ID	7378
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human UPP1 (Q16831).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16831
Protein Weight	34kDa
Shipping	Ice bag