


Anti-Uqcrfs1 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P47985

Item Number: CM0005813



Introduction

Recombinant rabbit monoclonal antibody against human UQCRC1, validated for WB, IF/ICC, ELISA. Detects human, mouse, rat. Suitable for studying mitochondrial complex III and electron transport chain.

[Learn More](#)

Parameter	Value
Product Name	UQCRC1 Rabbit mAb
Remarks/Alias	RIP1; RIS1; RISP; UQCRC5; MC3DN10; UQCRC1
Species	Human
GeneID (Human)	7386
GeneID	7386
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 175-274 of human UQCRC1 (P47985).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P47985
Protein Weight	30kDa
Shipping	Ice bag