


Anti-Usp16 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y5T5

Item Number: CM0029938



Introduction

Rabbit polyclonal antibody targeting human USP16 (94 kDa), validated for WB and ELISA. Cross-reacts with mouse, rat, and monkey. Suitable for studying histone H2A deubiquitination, mitosis, and gene regulation.

[Learn More](#)

Parameter	Value
Product Name	USP16 Rabbit pAb
Remarks/Alias	UBPM; UBP-M; USP16
Species	Human
GeneID (human)	10600
GeneID	10600
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human USP16 (NP_001001992.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat, Monkey
SWISS	Q9Y5T5
Protein Weight	94kDa
Shipping	Ice bag