

Anti-B-Amyloid(1-42) Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Db, Elisa - P05067

Item Number: CM0004994



Introduction

Rabbit monoclonal antibody against human β -Amyloid(1-42) for WB, IF-P, IHC-P, DB, ELISA. Detects a 4 kDa peptide, cross-reacts with mouse. Ideal for Alzheimer's research, diagnostic development, and mechanistic studies.

[Learn More](#)

Parameter	Value
Product Name	β -Amyloid(1-42) Rabbit mAb
Remarks/Alias	AAA; AD1; PN2; ABPP; APPI; CVAP; ABETA; PN-II; preA4; CTFgamma; alpha-sAPP; β -Amyloid(1-42)
Species	Human
Gene ID (Human)	351
Gene ID	351
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 691-770 of human β -Amyloid (NP_000475.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, DB, ELISA
Cross Reactivity	Mouse
SWISS	P05067
Protein Weight	4 kDa
Shipping	Ice bag