

Anti-B-Arrestin1 Rabbit Monoclonal Antibody For Wb And Elisa - P49407

Item Number: CM0003072



Introduction

High-quality rabbit monoclonal antibody against human β -arrestin1 (ARRB1) for Western blot and ELISA. Validated in human and mouse samples. Ideal for GPCR signaling research.

[Learn More](#)

Parameter	Value
Product Name	β -arrestin1 Rabbit mAb
Remarks/Alias	ARB1; ARR1; β -arrestin1
Species	Human
GeneID (Human)	408
GeneID	408
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human β -arrestin1 (P49407)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P49407
Protein Weight	47kDa
Shipping	Ice bag