


Anti-B1-Adrenergic Receptor Rabbit Polyclonal Antibody For Wb, Elisa - P08588

Item Number: CM0022732



Introduction

High-quality rabbit polyclonal antibody against human β 1-adrenergic receptor (ADRB1), validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for cardiovascular and neurological research in GPCR signaling and cardiac function studies.

[Learn More](#)

Parameter	Value
Product Name	β 1-Adrenergic Receptor Rabbit pAb
Remarks/Alias	RHR; B1AR; FNSS2; ADRB1R; BETA1AR; β 1-Adrenergic Receptor
Species	Human
Gene ID (Human)	153
Gene ID	153
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 350-450 of human β 1-Adrenergic Receptor (NP_000675.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08588
Protein Weight	51kDa
Shipping	Ice bag