

Anti-Vamp3/Cellubrevin Polyclonal Antibody For Wb And Elisa - Q15836

Item Number: CM0026881

Introduction

Rabbit polyclonal antibody against VAMP3/Cellubrevin, suitable for WB and ELISA. Recognizes human, mouse, and rat VAMP3. Targets the N-terminal domain and detects an 11kDa band.

[Learn More](#)

Parameter	Value
Product Name	VAMP3/Cellubrevin Rabbit pAb
Remarks/Alias	CEB; VAMP3/Cellubrevin
Species	Human
Gene ID (Human)	9341
Gene ID	9341
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-80 of human VAMP3/Cellubrevin (NP_004772.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15836
Protein Weight	11kDa
Shipping	Ice bag