


Anti-Vapb Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O95292

Item Number: CM0029769



Introduction

VAPB rabbit polyclonal antibody suitable for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat VAPB. For studying ER-endosome contact sites and UPR.

[Learn More](#)

Parameter	Value
Product Name	VAPB Rabbit pAb
Remarks/Alias	ALS8; VAP-B; VAMP-B; VAPB
Species	Human
GeneID	9217
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 10-222 of human VAPB (NP_004729.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95292
Protein Weight	27kDa
Shipping	Ice bag