


Anti-Vav2 Rabbit Polyclonal Antibody For Wb And Elisa - P52735

Item Number: CM0022333



Introduction

Anti-VAV2 rabbit polyclonal antibody validated for WB and ELISA, cross-reacting with human, mouse, and rat. VAV2 is a Rho GEF involved in angiogenesis.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] VAV2 Rabbit pAb
Remarks/Alias	VAV-2; V2
Species	Human
GeneID (human)	7410
GeneID	7410
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 570-700 of human VAV2 (NP_001127870.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P52735
Protein Weight	101kDa
Shipping	Ice bag