

Anti-Vdac1 Polyclonal Antibody (Kd Validated) For Wb, Elisa - P21796

Item Number: CM0024568



Introduction

KD validated rabbit polyclonal antibody targeting VDAC1 (voltage-dependent anion channel 1), suitable for WB and ELISA with human, mouse, rat reactivity. Essential for studying mitochondrial transport and apoptosis.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] VDAC1 Rabbit pAb
Remarks/Alias	PORIN; VDAC-1; [KD Validated] VDAC1
Species	Human
GeneID (human)	7416
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-283 of human VDAC1 (NP_003365.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P21796
Protein Weight	31kDa
Shipping	Ice bag