

Anti-Vdac1 / Porin Polyclonal Antibody For Wb, Ihc-P, Elisa - P21796

Item Number: CM0012835



Introduction

Rabbit polyclonal antibody against human VDAC1/Porin, validated for Western blot, IHC-P, and ELISA. Reacts with human, mouse, and rat samples. Ideal for mitochondrial and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	VDAC1 / Porin Rabbit pAb
Remarks/Alias	PORIN; VDAC-1; VDAC1 / Porin
Species	Human
Gene ID (Human)	7416
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-283 of human VDAC1 / Porin (NP_003365.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P21796
Protein Weight	31 kDa
Shipping	Ice bag