

Anti-Vegfa Polyclonal Antibody For Wb, If/Icc, Elisa - P15692

Item Number: CM0012979



Introduction

High-quality rabbit polyclonal antibody against VEGFA (VPF, VEGF) validated for WB, IF/ICC, ELISA in human, mouse, rat. Ideal for angiogenesis and hypoxia research. - P15692

[Learn More](#)

Parameter	Value
Product Name	VEGFA Rabbit pAb
Remarks/Alias	VPF; VEGF; MVCD1; VEGFA
Species	Human
GeneID (Human)	7422
GeneID	7422
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 91-191 of human VEGFA (NP_001165097.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15692
Protein Weight	27kDa
Shipping	Ice bag