


Anti-Vegfr1/Flt-1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P17948

Item Number: CM0014324



Introduction

High-quality rabbit polyclonal antibody targeting VEGFR1 (FLT-1). Validated for Western blot, immunofluorescence, and ELISA applications. Shows cross-reactivity with human, mouse, and rat samples. Ideal for angiogenesis, vasculature development, and cancer research studies.

[Learn More](#)

Parameter	Value
Product Name	VEGFR1 Rabbit pAb
Remarks/Alias	FLT; FLT-1; VEGFR1; VEGFR-1
Species	Human
GeneID (human)	2321
GeneID	2321
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1100-1180 of human VEGFR1 (NP_002010.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17948
Protein Weight	151kDa
Shipping	Ice bag