

Anti-Vglut2 Rabbit Monoclonal Antibody For Wb, If/Icc, Ihc-P, And Elisa - Q9P2U8

Item Number: CM0009910

Introduction

Recombinant rabbit monoclonal antibody against human VGLUT2. Validated for WB, IF/ICC, IHC-P, and ELISA in human, mouse, and rat tissues. Ideal for studying glutamate transport in synaptic vesicles and excitatory neurons.

[Learn More](#)

Parameter	Value
Product Name	VGLUT2 Rabbit mAb
Remarks/Alias	DNPI; VGLUT2
Species	Human
GeneID (Human)	57084
GeneID	57084
Immunogen	Synthetic peptide corresponding to amino acids 483-582 of human VGLUT2 (NP_065079.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-F, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9P2U8
Protein Weight	64kDa
Shipping	Ice bag