


Anti-Vhl Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P40337

Item Number: CM0022103



Introduction

Rabbit polyclonal antibody against human VHL tumor suppressor protein. Validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for ubiquitination and HIF pathway studies.

[Learn More](#)

Parameter	Value
Product Name	VHL Rabbit pAb
Remarks/Alias	RCA1; VHL1; pVHL; HRCA1; VHL
Species	Human
Gene ID (Human)	7428
Gene ID	7428
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human VHL (NP_937799.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P40337
Protein Weight	24kDa
Shipping	Ice bag