


Anti-Vps37A Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Nez2

Item Number: CM0025488



Introduction

VPS37A Rabbit Polyclonal Antibody validated for Western blot, IF/ICC, and ELISA. Reacts with human, mouse, and rat VPS37A, an ESCRT-I complex subunit involved in endosomal sorting and linked to liver cancer. Ideal for trafficking and disease research.

[Learn More](#)

Parameter	Description
Product Name	VPS37A Rabbit pAb
Remarks/Alias	HCRP1; PQBP2; SPG53; VPS37A
Species	Human
Gene ID (Human)	137492
Gene ID	137492
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 224-397 of human VPS37A (NP_689628.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NEZ2
Protein Weight	44kDa
Shipping	Ice bag