

Anti-Vti1A Rabbit Polyclonal Antibody For Wb, Elisa - Q96Aj9

Item Number: CM0029417

Introduction

High-quality rabbit polyclonal antibody against human VTI1A for WB and ELISA. Cross-reacts with mouse and rat. Suitable for vesicle transport, SNARE, and Golgi trafficking research. - Q96AJ9

[Learn More](#)

Parameter	Value
Product Name	VTI1A Rabbit pAb
Remarks/Alias	MMD53; MVti1; VTI1RP2; Vti1-rp2; VTI1A
Species	Human
GeneID (Human)	143187
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human VTI1A (NP_660207.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96AJ9
Protein Weight	25kDa
Shipping	Ice bag