


Anti-Wdr4 Rabbit Monoclonal Antibody For Wb, Ihc-P, Chip, Elisa - P57081

Item Number: CM0003391



Introduction

Recombinant rabbit monoclonal antibody against human WDR4 (P57081), validated for WB, IHC-P, ChIP, and ELISA. Targets the non-catalytic subunit of the METTL1 methyltransferase complex involved in tRNA/mRNA modification. Cross-reacts with mouse and rat.

[Learn More](#)

Parameter	Value
Product Name	WDR4 Rabbit mAb
Remarks/Alias	hWH, Wuho, MIGSB, TRM82, GAMOS6, TRMT82, WDR4
Species	Human
Gene ID	10785
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 250-350 of human WDR4 (P57081).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ChIP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P57081
Protein Weight	45kDa
Shipping	Ice bag