

Anti-Wdr45 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y484

Item Number: CM0026363



Introduction

Rabbit polyclonal antibody targeting WDR45, a key autophagy protein involved in autophagosome formation. Suitable for WB, IF/ICC, and ELISA in human, mouse, and rat samples. Ideal for autophagy research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	WDR45 Rabbit pAb
Remarks/Alias	JM5; NBIA4; NBIA5; WDRX1; WIPI4; WIPI-4; WDR45
Species	Human
Gene ID (Human)	11152
Gene ID	11152
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-148 of human WDR45 (NP_001025067.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9Y484
Protein Weight	40kDa
Shipping	Ice bag