

Anti-Wisp2 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - 076076

Item Number: CM0026231

Introduction

Rabbit anti-WISP2 polyclonal antibody validated for WB, IHC-P, IF/ICC, ELISA. Targets human WISP2 (CCN5), a secreted protein involved in bone turnover and cell adhesion. Cross-reacts with mouse and rat.

[Learn More](#)

Parameter	Value
Product Name	WISP2 Rabbit pAb
Remarks/Alias	CT58; WISP2; CTGF-L
Species	Human
GeneID	8839
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 164-250 of human WISP2 (NP_003872.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	076076
Protein Weight	27kDa
Shipping	Ice bag