

Anti-Wtap Monoclonal Antibody For Wb, Ihc-P, Elisa - Q15007

Item Number: CM0004348



Introduction

WTAP Rabbit monoclonal antibody for WB, IHC-P, ELISA. Recognizes human, mouse, rat. Targets pre-mRNA-splicing regulator involved in m6A methylation and cell cycle. Ideal for IVD research.

[Learn More](#)

Parameter	Specification
Product Name	WTAP Rabbit mAb
Remarks/Alias	Mum2; WTAP
Species	Human
GeneID	9589
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 212-355 of human WTAP (NP_004897.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15007
Protein Weight	44kDa
Shipping	Ice bag