


## Anti-Xpnpep1 Polyclonal Antibody For Wb, Elisa - Q9Nqw7

Item Number: CM0013806



### Introduction

XPNPEP1 Rabbit pAb is a polyclonal antibody targeting human Xaa-Pro aminopeptidase 1, validated for WB and ELISA. It cross-reacts with mouse, and detects a 70 kDa cytosolic protein involved in peptide degradation and bradykinin metabolism.

[Learn More](#)

| Parameter        | Value                                                                                                                |
|------------------|----------------------------------------------------------------------------------------------------------------------|
| Product Name     | XPNPEP1 Rabbit pAb                                                                                                   |
| Remarks/Alias    | APP1; SAMP; XPNPEP; XPNPEPL; XPNPEPL1; XPNPEP1                                                                       |
| Species          | Human                                                                                                                |
| Gene ID (Human)  | 7511                                                                                                                 |
| Gene ID          | 7511                                                                                                                 |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 44-240 of human XPNPEP1 (NP_065116.3). |
| Source           | Rabbit                                                                                                               |
| Category         | Polyclonal Antibodies                                                                                                |
| Application      | WB, ELISA                                                                                                            |
| Cross Reactivity | Human, Mouse                                                                                                         |
| SWISS            | Q9NQW7                                                                                                               |
| Protein Weight   | 70kDa                                                                                                                |
| Shipping         | Ice bag                                                                                                              |