


## Anti-Xrcc1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P18887

Item Number: CM0013812



### Introduction

Rabbit polyclonal antibody against XRCC1, validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Suitable for DNA repair and single-strand break research.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | XRCC1 Rabbit pAb                                                                                                   |
| Remarks/Alias    | RCC; SCAR26; XRCC1                                                                                                 |
| Species          | Human                                                                                                              |
| GeneID (human)   | 7515                                                                                                               |
| GeneID           | 7515                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 480-633 of human XRCC1 (NP_006288.2) |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                                           |
| Cross Reactivity | Human, Mouse, Rat                                                                                                  |
| SWISS            | P18887                                                                                                             |
| Protein Weight   | 69kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |