


Anti-Ykl-40/Chi3L1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q61362

Item Number: CM0028504



Introduction

Rabbit polyclonal antibody against mouse YKL-40/CHI3L1 (CHI3L1; BRP39; GP39). Validated for WB, IHC-P, ELISA. Cross-reacts with human. Involved in inflammation, tissue remodeling, and antimicrobial responses.

[Learn More](#)

Parameter	Value
Product Name	YKL-40/CHI3L1 Rabbit pAb
Remarks/Alias	GP39; ASRT7; GP-39; YK-40; YKL40; CGP-39; YKL-40; YYL-40; HC-gp39; HCGP-3P; hCGP-39
Species	Mouse
GeneID (human)	1116
GeneID	12654
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-389 of mouse YKL-40/CHI3L1 (NP_031721.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q61362
Protein Weight	43 kDa
Shipping	Ice bag