

Anti-Yme1L1 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q96Ta2

Item Number: CM0009890

Introduction

YME1L1 Rabbit Monoclonal Antibody validated for WB, IF/ICC, and ELISA. Detects human YME1L1, an ATP-dependent zinc metalloprotease that regulates mitochondrial morphology, protein quality control, and stress adaptation.

[Learn More](#)

Parameter	Value
Product Name	YME1L1 Rabbit mAb
Remarks/Alias	FTSH; MEG4; PAMP; OPA11; YME1L
Species	Human
Gene ID (Human)	10730
Gene ID	10730
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 580-773 of human YME1L1 (NP_647473.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q96TA2
Protein Weight	86kDa
Shipping	Ice bag