

## Anti-Yy2 Polyclonal Antibody For Wb, If/Icc, Elisa - O15391

Item Number: CM0023874

### Introduction

Rabbit polyclonal antibody against YY2 (ZNF631), validated for WB, IF/ICC, and ELISA in human, mouse, rat. Ideal for transcription regulation and development research.

[Learn More](#)

| Parameter        | Value                                                                                                                    |
|------------------|--------------------------------------------------------------------------------------------------------------------------|
| Product Name     | YY2 Rabbit pAb                                                                                                           |
| Remarks/Alias    | ZNF631; YY2                                                                                                              |
| Species          | Human                                                                                                                    |
| GeneID (human)   | 404281                                                                                                                   |
| GeneID           | 404281                                                                                                                   |
| Immunogen        | Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 200-300 of human YY2 (NP_996806.2). |
| Source           | Rabbit                                                                                                                   |
| Category         | Polyclonal Antibodies                                                                                                    |
| Application      | WB, IF/ICC, ELISA                                                                                                        |
| Cross Reactivity | Human, Mouse, Rat                                                                                                        |
| SWISS            | O15391                                                                                                                   |
| Protein Weight   | 41kDa                                                                                                                    |
| Shipping         | Ice bag                                                                                                                  |