


Anti-Zbtb48 Polyclonal Antibody For Wb, Elisa - P10074

Item Number: CM0029129



Introduction

Rabbit polyclonal antibody against human ZBTB48 (TZAP), suitable for WB and ELISA, cross-reacts with mouse. Useful for investigating telomere length regulation and gene expression.

[Learn More](#)

Parameter	Value
Product Name	ZBTB48 Rabbit pAb
Remarks/Alias	HKR3; TZAP; ZNF855; pp9964; ZBTB48
Species	Human
GeneID (Human)	3104
GeneID	3104
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human ZBTB48 (NP_005332.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P10074
Protein Weight	77kDa
Shipping	Ice bag