


Anti-Zeb1 Polyclonal Antibody For Wb, Ihc-P, Elisa - P37275

Item Number: CM0014416



Introduction

Rabbit polyclonal antibody to ZEB1, a zinc finger transcription factor involved in epithelial-mesenchymal transition. Validated for WB, IHC-P, and ELISA in human, mouse, and rat samples. Molecular weight: 124 kDa.

[Learn More](#)

Parameter	Specification
Product Name	ZEB1 Rabbit pAb
Remarks/Alias	BZP; TCF8; AREB6; FECD6; NIL2A; PPCD3; ZFHEP; ZFH1A; DELTAEF1; ZEB1
Species	Human
Gene ID (Human)	6935
Gene ID	6935
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 490-710 of human ZEB1 (NP_110378.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P37275
Protein Weight	124kDa
Shipping	Ice bag