


Anti-Zmat2 Rabbit Polyclonal Antibody For Wb, Elisa - Q96Nc0

Item Number: CM0015341



Introduction

Rabbit polyclonal antibody targeting ZMAT2, a spliceosome component involved in pre-mRNA splicing. Validated for WB and ELISA applications in human, mouse, and rat samples. Ideal for splicing and gene expression research.

[Learn More](#)

Parameter	Value
Product Name	ZMAT2 Rabbit pAb
Remarks/Alias	Snu23; Ptg-12; ZMAT2
Species	Human
Gene ID (Human)	153527
Gene ID	153527
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-199 of human ZMAT2 (NP_653324.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96NC0
Protein Weight	24kDa
Shipping	Ice bag