

Anti-Zmsnf2 Rabbit Polyclonal Antibody For Wb, Elisa - A0A3L6Fqf6

Item Number: CM0018771

Introduction

Rabbit polyclonal antibody targeting maize ZmSNF2 (UniProt A0A3L6FQF6), validated in WB and ELISA. Ideal for chromatin remodeling research in maize. Recombinant protein immunogen, shipped with ice bag.

[Learn More](#)

Parameter	Value
Product Name	ZmSNF2 Rabbit pAb
Remarks/Alias	GRMZM2G167668; GRMZM2G379804; ZmSNF2
Species	Maize
Gene ID	103652266
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 901-1089 of maize ZmSNF2 XP_020405944.1
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Maize
SWISS	A0A3L6FQF6
Protein Weight	155kDa
Shipping	Ice bag