

## Anti-Znf20 Monoclonal Antibody For Wb, If/Icc, If-P, Elisa - P17024

Item Number: CM0008078

### Introduction

Anti-ZNF20 rabbit monoclonal antibody validated for WB, IF/ICC, IF-P, ELISA. Recognizes human, mouse, rat. Target is ZNF20 (KOX13), a 62kDa zinc finger transcription factor. Ideal for gene regulation research.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | ZNF20 Rabbit mAb                                                                                                  |
| Remarks/Alias    | KOX13                                                                                                             |
| Species          | Human                                                                                                             |
| GeneID (Human)   | 7568                                                                                                              |
| GeneID           | 7568                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 51-118 of human ZNF20 (NP_066966.2) |
| Source           | Rabbit                                                                                                            |
| Category         | Monoclonal Antibodies                                                                                             |
| Application      | WB, IF/ICC, IF-P, ELISA                                                                                           |
| Cross Reactivity | Human, Mouse, Rat                                                                                                 |
| SWISS            | P17024                                                                                                            |
| Protein Weight   | 62kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                           |