

Anti-Znf331 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Nqx6

Item Number: CM0015096

Introduction

Anti-ZNF331 rabbit polyclonal antibody validated for WB and ELISA; detects ZNF331 in human, mouse, and rat. Targets a testis-specific transcription factor potentially involved in spermatogenesis. UniProt ID: Q9NQX6.

[Learn More](#)

Parameter	Value
Product Name	ZNF331 Rabbit pAb
Remarks/Alias	RITA; ZNF361; ZNF463; ZNF331
Species	Human
Gene ID (Human)	55422
Gene ID	55422
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-130 of human ZNF331 (NP_061025.5).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NQX6
Protein Weight	54kDa
Shipping	Ice bag