

Anti-Znf521 Polyclonal Antibody For Wb, Elisa - Q96K83

Item Number: CM0021004

Introduction

Rabbit polyclonal antibody targeting human ZNF521 (EHZF/Evi3) for Western blot and ELISA applications. Demonstrates cross-reactivity with human and rat samples. Ideal for stem cell and transcriptional regulation research.

[Learn More](#)

Parameter	Value
Product Name	ZNF521 Rabbit pAb
Remarks/Alias	EHZF; Evi3; ZNF521
Species	Human
GeneID (Human)	25925
GeneID	25925
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 210-420 of human ZNF521 (NP_056276.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q96K83
Protein Weight	148kDa
Shipping	Ice bag