

Anti-Znf785 Polyclonal Antibody For Wb, Elisa - A8K8V0

Item Number: CM0030861



Introduction

High-quality rabbit polyclonal antibody against human ZNF785, reactive with human and mouse, for WB and ELISA applications. Validated for research use in transcriptional regulation studies.

[Learn More](#)

Parameter	Value
Product Name	ZNF785 Rabbit pAb
Remarks/Alias	ZNF688; ZNF785
Species	Human
Gene ID (Human)	146540
Gene ID	146540
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 200-405 of human ZNF785 (NP_689671.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	A8K8V0
Protein Weight	46kDa
Shipping	Ice bag