


Anti-Zo-1 Rabbit Monoclonal Antibody For Wb, If/Icc, If-F, Ihc-P, Elisa - P39447

Item Number: CM0005134



Introduction

Recombinant rabbit monoclonal antibody targeting mouse ZO-1 tight junction protein. Validated for WB, IF/ICC, IF-F, IHC-P, ELISA; cross-reacts with human, mouse, rat. ~195 kDa. Suitable for studying tight junction biology.

[Learn More](#)

Parameter	Value
Product Name	ZO-1 Rabbit mAb
Remarks/Alias	ZO1
Species	Mouse
GeneID (Human)	7082
GeneID (Mouse)	21872
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1525-1621 of mouse ZO-1 (NP_033412.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-F, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P39447
Protein Weight	195 kDa
Shipping	Ice bag