

Anti-Zp3 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P21754

Item Number: CM0016425

Introduction

High-quality ZP3 Rabbit pAb for WB, IF-P, ELISA applications. Recognizes human, mouse, and rat ZP3 protein, a key sperm receptor in the zona pellucida. Ideal for fertility and reproductive biology studies.

[Learn More](#)

Parameter	Value
Product Name	ZP3 Rabbit pAb
Remarks/Alias	ZPC; ZP3A; ZP3B; Zp-3; OOMD3; OZEMA3; ZP3
Species	Human
GeneID (Human)	7784
GeneID	7784
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 246-345 of human ZP3 (NP_001103824.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P21754
Protein Weight	47kDa
Shipping	Ice bag