

Anti-Zyxin Monoclonal Antibody For Wb, If/Icc, Elisa - Q15942

Item Number: CM0006641

Introduction

Recombinant rabbit monoclonal antibody against human Zyxin (Q15942), suitable for WB, IF/ICC, ELISA. Detects endogenous protein at 61kDa. Cross-reacts with mouse. Ideal for cell adhesion and focal adhesion research.

[Learn More](#)

Parameter	Value
Product Name	Zyxin Rabbit mAb
Remarks/Alias	ESP-2; HED-2; Zyxin
Species	Human
GeneID (Human)	7791
GeneID	7791
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 250-350 of human Zyxin (Q15942).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15942
Protein Weight	61kDa
Shipping	Ice bag