


Antxr1 Rabbit Pab For Wb And Elisa - Q9H6X2

Item Number: CM0029550



Introduction

ANTXR1 Rabbit Polyclonal Antibody validated for WB, ELISA; cross-reacts with Human, Mouse. Targets anthrax toxin receptor 1/TEM8, a cell adhesion and migration protein. Ideal for IVD research.

[Learn More](#)

Parameter	Value
Product Name	ANTXR1 Rabbit pAb
Remarks/Alias	ATR; GAPO; TEM8; ANTXR1
Species	Human
GeneID (Human)	84168
GeneID	84168
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 84-333 of human ANTXR1 (NP_060623.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9H6X2
Protein Weight	63kDa
Shipping	Ice bag