


Aoc1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P19801

Item Number: CM0029314



Introduction

Rabbit polyclonal antibody raised against human AOC1 (diamine oxidase/histaminase). Verified for WB, IF/ICC, ELISA; cross-reacts with mouse and rat. Ideal for allergy, immune response, and amine metabolism research.

[Learn More](#)

Parameter	Value
Product Name	AOC1 Rabbit pAb
Remarks/Alias	ABP; DAO; KAO; ABP1; DAO1; AOC1
Species	Human
GeneID	26 (human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-300 of human AOC1 (NP_001082.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P19801
Protein Weight	85kDa
Shipping	Ice bag