


Aromatase (Cyp19A1) Rabbit Pab For Wb, Ihc-P, If/Icc, Elisa - P11511

Item Number: CM0015257



Introduction

Validated for WB, IHC-P, IF/ICC, and ELISA, this rabbit polyclonal antibody targets human Aromatase (CYP19A1), a key enzyme in estrogen biosynthesis. Cross-reacts with human and rat. Ideal for hormone-dependent cancer and reproductive biology research.

[Learn More](#)

Parameter	Value
Product Name	Aromatase (CYP19A1) Rabbit pAb
Remarks/Alias	ARO; ARO1; CPV1; CYAR; CYP19; CYPXIX; P-450AROM; Aromatase (CYP19A1)
Species	Human
GeneID (Human)	1588
GeneID	1588
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human Aromatase (CYP19A1) (NP_000094.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	P11511
Protein Weight	58kDa
Shipping	Ice bag