

Asymmetric Dimethyl-Histone H3-R8 Rabbit Pab (Q16695 / P68431) - Validated For Wb, Ihc-P, Chip, Elisa

Item Number: CM0000374



Introduction

Asymmetric DiMethyl-Histone H3-R8 rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, ELISA, and ChIP. Recognizes human, mouse, rat, and other species (wide range predicted). Ideal for epigenetic studies of arginine methylation and chromatin regulation. UniProt: Q16695 / P68431.

[Learn More](#)

Parameter	Value
Product Name	Asymmetric DiMethyl-Histone H3-R8 Rabbit pAb
Remarks/Alias	H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3; Asymmetric DiMethyl-Histone H3-R8
Species	Human
GeneID (Human)	H3
GeneID	8290/8350
Immunogen	Synthetic peptide - A synthetic asymmetric dimethylated peptide around R8 of human histone H3 (NP_003520.1)
Source	Rabbit
Category	Methylated Antibodies
Application	WB, IHC-P, IF/ICC, ELISA, ChIP
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag