


## Cadm1 Rabbit Pab For Wb, Elisa - Q9By67

Item Number: CM0023274



### Introduction

CADM1 Rabbit pAb is a polyclonal antibody validated for WB and ELISA. Reacts with human, mouse, rat CADM1. Ideal for studying cell adhesion, tumor suppression, and synaptic function.

[Learn More](#)

| Parameter        | Value                                                                                                                                   |
|------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | CADM1 Rabbit pAb                                                                                                                        |
| Remarks/Alias    | BL2; ST17; IGSF4; NECL2; RA175; TSLC1; IGSF4A; Necl-2; SYNCAM; sgIGSF; sTSLC-1; synCAM1; CADM1                                          |
| Species          | Human                                                                                                                                   |
| Gene ID (Human)  | 23705                                                                                                                                   |
| Gene ID          | 23705                                                                                                                                   |
| Immunogen        | Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 45-344 of human CADM1 (NP_055148.3). |
| Source           | Rabbit                                                                                                                                  |
| Category         | Polyclonal Antibodies                                                                                                                   |
| Application      | WB, ELISA                                                                                                                               |
| Cross Reactivity | Human, Mouse, Rat                                                                                                                       |
| SWISS            | Q9BY67                                                                                                                                  |
| Protein Weight   | 49kDa                                                                                                                                   |
| Shipping         | Ice bag                                                                                                                                 |