

Calmegin Rabbit Mab For Wb, If-P, Elisa - O14967

Item Number: CM0003867

Introduction

Calmegin Rabbit Monoclonal Antibody validated for WB, IF-P, ELISA. Recognizes human, mouse, rat Calmegin (70kDa). For spermatogenesis and ER chaperone studies. - O14967

[Learn More](#)

Parameter	Value
Product Name	Calmegin Rabbit mAb
Species	Human
GeneID (Human)	1047
GeneID	1047
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 180-350 of human Calmegin (O14967).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14967
Protein Weight	70kDa
Shipping	Ice bag