


Capg Rabbit Monoclonal Antibody For Wb, Elisa - P40121

Item Number: CM0009574



Introduction

CAPG Rabbit Monoclonal Antibody (P40121) validated for Western blot and ELISA. Reacts with human, mouse, rat. Targets macrophage-capping protein, a key actin regulator. Suitable for research on cytoskeleton and macrophage function.

[Learn More](#)

Parameter	Value
Product Name	CAPG Rabbit mAb
Remarks/Alias	MCP; AFCP; HEL-S-66; CAPG
Species	Human
GeneID (Human)	822
GeneID	822
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-188 of human CAPG (P40121).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P40121
Protein Weight	38kDa
Shipping	Ice bag