

Cbfb Rabbit Pab For Wb, Ihc-P, If/Icc, Elisa - Q13951

Item Number: CM0028808



Introduction

CBFB rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, ELISA applications. Recognizes human, mouse, and rat core-binding factor subunit beta, a key transcriptional coregulator in hematopoiesis and immune function. Suitable for cancer and infectious disease research.

[Learn More](#)

Parameter	Value
Product Name	CBFB Rabbit pAb
Remarks/Alias	CLCD2; PEBP2B; CBFB
Species	Human
GeneID (Human)	865
GeneID	865
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-182 of human CBFB (NP_001746.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13951
Protein Weight	22kDa
Shipping	Ice bag