

Cd200 Rabbit Pab For Wb, Elisa - P41217

Item Number: CM0026524

Introduction

High-quality CD200 Rabbit pAb (Polyclonal Antibodies) for WB, ELISA. Targets human CD200 (OX-2 membrane glycoprotein; SWISS: P41217), cross-reacts with mouse and rat. Suitable for investigating T-cell costimulation and myeloid cell regulation.

[Learn More](#)

Parameter	Value
Product Name	CD200 Rabbit pAb
Remarks/Alias	MRC; MOX1; MOX2; OX-2; CD200
Species	Human
Gene ID (Human)	4345
Gene ID	4345
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-240 of human CD200 (NP_005935.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P41217
Protein Weight	31kDa
Shipping	Ice bag